

Technical Data Sheet

Doc No – Ab/TDS/LH-65

Rev No : 03

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Rava Pipette Controller, PP

Valid for Catalog No. E11113

1. Description/Specification

1.1 Description:

Rave Pipette controller ideal for Repetitive Pipetting performance

Maximizing Pipetting comfort

1.2 Design:

Robust, Light weight & simple to use operation.

1.3 Performance Specifications:

1.3.1 Can use volumetric or serological pipettes ranging from 0.1 ml, to 100 ml.

1.3.2 It contains integral & replacement 0.3 μ m Filter.

1.3.3 Autoclavable silicon pipette holder

1.4 Material : PP

1.5 Color : Grey

1.6 Basic Fetures

1.6.1 Made up of highly Virgin thermoplastic polypropylene

1.6.2 Autoclavability : Yes (Only Silicon Pipette Holder – Don't autoclave whole instrument)



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos


Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com

Technical Data Sheet

1.7 Dimension :

Code	Dimension
E11113	: Length- 210 mm, ϕ -60mm
E11119	: (Spare) 0.3 μ m Filter



Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata - 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com