

Technical Data Sheet

Doc No – Ab/TDS/RB-37

Rev No : 03

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Rotator Work Station, PP

Valid for Catalog No. P20208

1. Description/Specification

1.1 Description:

Reusable PP (Poly Propylene) Rotator Work Station

1.2 Design:

Cube Type Multi Hole for keeping Different Product as Per own requirement

1.3 Performance Specification:

1.3.1 Product	No. per unit
50 ml Centrifuge Tubes	2 per unit
15 ml Centrifuge Tube	8 per unit
1.5 & 2.5 Micro Centrifuge Tube	9 Per unit
0.5 ml Micro centrifuge tube	6 per unit
PCR Plates (Unskirted & semi Skirted 96 well)	1 per 2 units put together
PCR 8 Tubes strip	6 per unit
PCR 12 Tubes Strip	8 per 2 unit
PCR 0.2 ml Tubes	48 per unit

Universal fit all above Product as mentioned.

1.4 Material:

1.4.1 Made up with high impact Strength poly propylene

1.4.2 Compile with requirement of directive 67/548/EEC or 99/45/EC

1.4.3 Compile with requirement of Directive 2002/72/EC and requirement of FDA Regulation 21 CFR 177.1520. olefin polymers.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

1.4.4 Complies with the requirement of Indian Standard IS: 10910 on “ Specification for Polypropylene and its copolymer for its safe use in contact with food stuffs, Pharmaceutical and Drinking water”

1.5 Color :

Assorted Colour - 4 Pcs per case – Blue, Green, Yellow, Red

1.6 Sterilization :

Can be autoclaved at 15 psi, 121⁰C for 20 minutes

1.7 Quality control:

Testing of various characteristic in accordance with the valid Specification

1.8. Basic Features:

1.8.1 Made up of High pure virgin Polypropylene (PP)

1.8.2 Temperature range: -80⁰C to +121⁰C

1.8.3 Autoclavability: Yes, Autoclavable at 15 psi and 121⁰C for 20 minutes

1.8.4 Chemical resistance: See website.

1.8.5 Interlocking Racks – Convenient to place all type of product in single racking products

1.8.6 Customize Specific requirements according to end user.


1.8.7 Rotate each of the four modules independently


1.8.8 Rack maintains its orientation until twist again.

1.8.9 Dimension: 108 x 68 x108 mm

Autoclaving Guideline for Plastics

All plastic consumables are autoclavable at 121⁰C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50⁰C before removing product.


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com