

# Technical Data Sheet

Doc No – Ab/TDS/LH-72

Rev No : 02

## TECHNICAL DATA SHEET

**Abdos Labtech Pvt. Ltd.**

**Serological Pipettes Sterile (PS)**

**Valid for Catalog No. P10601, P10602, P10603, P10604, P10605, P10606**

### **1. Description/Specification**

#### **1.1 Description:**

Disposable Serological Pipettes Sterile (PS)

#### **1.2 Design:**

Transparent Cylindrical Design

#### **1.3 Performance Specifications:**

P10601 (1 ml), P10602 (2 ml), P10603 (5 ml), P10604 (10 ml), P10605 (25 ml), P10606 (50 ml)

Accuracy within  $\pm 2\%$  At Full volume.

#### **1.4 Material:**

1.4.1 Comply with requirement of Directive EU No 10/2011 and requirement of FDA Regulation 21 CFR 177.1640

1.4.2 Materials are manufactured according to the requirements of GMP Regulation (EC) No 2023/2006

1.4.3 Material does not contain any Hazardous chemical substance & compliance with REACH AVHC & RoHS Free.

#### **1.5 Color:**

Transparent (Natural)

#### **1.6 Graduations:**

Large Cylindrical with temperature, chemical and solvent resistant ink.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttarakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com

# Technical Data Sheet

## 1.7 Sterilization:

**Sterile catalog number** from this lot have been pre-sterilized and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products Requirements for validation and Routine control- Radiation sterilization*. The product has been sterilized with Co-60 Gamma Radiation

## 1.8 Quality control:

Testing of various characteristic in accordance with the valid Specification

## 1.9 other Information:

Recommended for single use only

## 2. Basic Features

**2.1 Temperature range:** -80<sup>0</sup>C to +90<sup>0</sup>C

**2.2 Autoclavability:** NO

**2.3 Chemical resistance:** See website.

## Autoclaving Guideline for Plastics

Polystyrene can't be autoclaved.



Shiv Ganga Industrial Estate,  
Plot no. 1 Lakeshwari, Bhagwanpur,  
Roorkee - 247667 Uttrakhand, India  
**An ISO 9001:2008 Company**  
A Division of Abdos

**Regd. Office :**  
13A, Government Place (East)  
Kolkata – 700069, India  
**CIN :** U85100WB2008PTC123903

**P :** +91 (0) 1332 235277  
**E :** info@abdoslabware.com  
**W :** www.abdoslabware.com