

Technical Data Sheet

Doc No–Ab/TDS/BC-98

Rev No: 00

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Wide Mouth Bottle with Sealing Cap, HDPE

Valid for Catalog No. P11201, P11202, P11203, P11204

1. Description/Specification

1.1 Description:

Wide Mouth Bottle with Sealing Cap, HDPE

1.2 Raw material:

1.2.1 Body: HDPE (High Density poly Ethylene)

1.2.2 Cap: HDPE (High Density poly Ethylene)

1.3 Other features:

Its use is recommended when you need to preserve the content after the filling. They are suitable for the storage of very fluid liquids. The complete emptying is guaranteed by the inclination of the shoulders.

1.4 DIMENSIONAL SPECIFICATIONS:

| Code | Capacity (ml) | Mouth Diameter (mm) | Height (mm) | No. per Case |
|---------------|----------------------|----------------------------|--------------------|---------------------|
| P11201 | 250 | 30 | 140 | 12 |
| P11202 | 500 | 30 | 179 | 6 |
| P11203 | 750 | 30 | 214 | 6 |
| P11204 | 1000 | 30 | 218 | 6 |

1.5 GENERAL FEATURES:



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com

Technical Data Sheet

| | |
|-----------------|------------------------------|
| CE mark | No |
| Autoclavability | No |
| Translucence | Yes |
| Resistant | To be effected by radiations |
| Sterilization | From – 40°C to + 80°C |

1.6 Chemical resistance: See website. **Autoclaving Guideline for Plastics.**

1.7 Other Information (not applicable if neglected):


1.7.1 fair resistance; exposure (30 days) to the chemical does not cause any damage.


1.7.2 sufficient resistance; exposure causes damage of poor importance, which sometimes is only temporary.

1.7.3 poor resistance; exposure to chemical is not allowed, or causes immediate permanent damage.

Results of testing with reference materials are to be considered as indication: in case of specific use it is recommended to carry out preliminary testing.

Information relating to the resistances of plastic materials to temperatures, sterilization and washing treatments, have been made as per reports of the raw materials' producers, following the data published in the literature and the experience gained in the use of products


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com