

Technical Data Sheet

Doc No– Ab/TDS/Cryo-30

Rev No: 01

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Freezing Cardboard Cryo Boxes, Cardboard

Valid for Catalog No. P20606, P20613

1. Description/Specification

1.1 Description: Freezing Cardboard Cryo Boxes, Cardboard

1.1.1 Freezing Cardboard Cryo Boxes designed for storage & preservation of frozen samples.

1.1.2 Freezing Cardboard Cryo Boxes can withstand ultra low temperature in deep freezers.

1.1.3 Freezing Cardboard Cryo Boxes are removable top lid makes it easy to open & close.

1.1.4 Lamination on these boxes does not allow moisture to soak or liquid to stick on these boxes.

1.2 Design:

Square Box with Lid

1.3 Performance Specification:

1.3.1 Capacity: P20606 – 1.5, 2.0 ml Micro tube and 1.0, 2.0 ml Cryo tubes – 81 Places (9x9 Array)

1.3.2 Capacity: P20613 – 1.5, 2.0 ml Micro tube and 1.0, 2.0 ml Cryo tubes – 100 Places (10x10 Array)

1.3.3 P20606 – Dimension (132x132x50)

1.3.4 P20613 - Dimension (132x132x50)

1.4 Color:

Box: White

1.5 Quality control:

Testing of various characteristic in accordance with the valid Specification



Technical Data Sheet

1.6 Other Information:

Reusable, can withstand ultra low temperature in deep freezers.

1.7 Basic Features

1.7.1 Made up of Cardboard.

1.7.2 Performance: Freezing Cardboard Boxes are removable top lid makes it easy to open & close & Lamination on these boxes does not allow moisture to soak or liquid to stick on these boxes.


1.7.3 Freezing Cardboard Boxes can withstand ultra low temperature in deep freezers.

1.8 Basic Features:

1.8.1 Temperature range: Tube:-196 °C to + 121°C

1.8.2 Autoclavability: Yes, Autoclavable at 15 psi and 121 °C for 20 minutes

1.8.3 Chemical resistance: see website.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com