

Technical Data Sheet

Doc No – Ab/TDS/CW-02

Rev No : 02

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Micro Pestles PP

Valid For Cat. No. P10206

1. Description/Specification

1.1 Description:

Disposable Micro Pestles (PP), DNase, RNase and Pyrogen Free

1.2 Design:

Tube: Cylindrical with Smooth conical Head Made up of Highest quality Virgin Polypropylene.

1.3 Performance Specifications:

P10206 Natural

1.4 Material:

1.4.1 **Tube:** Made up of High quality virgin polypropylene, confirm US FDA 21CFR

1.5 Color: Natural

1.6 Graduation and writing surface:

1.6.1 85 mm Tube with smooth surface

1.7 Sterilization:

DNase, RNase and Pyrogen Free

Sterile catalog number from this lot have been pre-sterilizes and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products Requirements for validation control- Radiation sterilization*. The product has been sterilized with Co-60 Gamma Radiation.

1.8 Quality control:

Testing of various characteristics in accordance with the valid Specification

1.9 Other information

Recommended for single use only

Pack Size: 10pieces/case



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

2. Basic Features


2.1 Temperature range: Tube:-80⁰C to +121⁰C

2.2 Autoclavability: Yes, Autoclavable at 15 psi and 121⁰C for 15 minutes

2.3 Chemical resistance: See website.

Autoclaving Guideline for Plastics

All plastic consumables are autoclavable at 121⁰C for 15 minutes, 15 psi/1atm. Do not use cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50⁰C before removing product.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com