

Technical Data Sheet

Doc No – Ab/TDS/Petri-91-I

Rev No: 01

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Petri Dish, PS - Aseptic

Valid for Cat. No. P10901, P10902, P10903, P10904

1. Description/Specification

1.1 Description: Petri Dish- Aseptic.

1.2 Design: Shallow dish design

1.3 Performance Specifications:

Code	Diameter (mm)	Height (mm)	Compartment	Vented	Packing	No. Per Case
P10901	90	15.80	No	3 Vent	Sleeve of 20	480
P10902	90	15.80	No	3 Vent	Individual Wrapped	380
P10903	90	15.80	No	0 Vent	Sleeve of 20	480
P10904	90	15.80	No	0 Vent	Individual Wrapped	380

1.4 Material:

1.4.1 Petri Dish made up of Hydrophobic Polystyrene (PS).

1.4.2 Confirm with US FDA Title 21CFR section '177.1640 of US Food and Drug Administration, regarding use of polystyrene in food contact application.

1.4.3 Meets with various stipulations laid down by European Union directive 20'11110/EU dated 1410112011 relating to Plastic materials and articles intended to come into contact with food and its amendments (Commission Regulation (EU) No. 1183 t2012).

1.4.4 Complies with the requirement of Indian Standard IS: 10910 on "Specification for Polypropylene and its copolymer for its safe use in contact with food stuffs, Pharmaceutical and Drinking water"

1.4.5 Reference to Clause (24) of 20111101EU, Dual use additives subject to restrictions in food are not intentionally used to produce this product.

1.5 Color: Transparent (Natural color)



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

1.6 Quality control:

DNase, RNase, and Pyrogen Free

Testing of various characteristics in accordance with the valid Specification

1.7 Other information

1.7.1 Petri Dishes are suitable for all bacteriological applications.

1.7.2 Vented dishes have elevated lid that allows for superior gas exchange and evaporation.


1.7.3 Petri dishes retain their shape up to 55°C.


1.8 Basic Features

1.8.1 Temperature range: Stable up to 55°C.

1.8.2 Autoclavability: NO

1.8.3 Chemical resistance: See website.


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata - 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com