

Technical Data Sheet

Doc No–Ab/TDS/UTL-160

Rev No: 00

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Safety Goggles, PC

Valid for Catalog No. U40101

1. Description/Specification

1.1 Description: Safety Goggles, PC

1.1.1 Polycarbonate lens with side shield.

1.1.2 Impact resistant and lightweight frames for comfortable eyewear.

1.1.3 99.9% UVA / UVB protection.

1.2 Design: Safety Goggles

1.3 Performance Specifications:

Code	No. Per Case
U40101	2

1.4 Material: PC (Polycarbonate)

1.5 Color: Transparent Natural

1.6 Quality control:


Testing of various characteristics in accordance with the valid Specification


1.7 Other information

1.7.1 Excellent clarity and visibility.

1.7.2 Moderate chemical stability.

1.7.3 Temperature range from -130°C to +125°C


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com