

Technical Data Sheet

Doc No – Ab/TDS/Eqpt-136

Rev No : 02

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

SWIRLMAX TUBE ROLLER

Valid for Catlog No. E11230

1. Description/Specification

1.1 Description:

SWIRLMAX Tube Roller can be used for rocking and rolling for complete mixing. Idle for mixing blood samples, viscous substances and liquid solid suspensions. Prevention of blood coagulation and latex diagnostics. It can be used in cold rooms and incubators.

1.2 Design :

Classis Type

1.3 performance Specifications :

- 1.3.1 Voltage (AC) : 220~240 Volt
- 1.3.2 Frequency (Hz) : 50
- 1.3.3 Power (W) : 25
- 1.3.4 Amplitude : 24 mm
- 1.3.5 No of rollers : 6
- 1.3.6 Roller size : 280mm
- 1.3.7 Motor type : DC brushles Motor
- 1.3.8 Speed range (RPM) : 0 ~ 70
- 1.3.9 Run type : Continuous operation
- 1.4.0 overall dimension (mm) : 450 x 260 x 120 mm
- 1.4.1 Weight : 4.5 Kg
- 1.4.2 Permissible ambient Temp : 5- 40
- 1.4.3 Permissible Relative humadity : 80 %



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos


Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903


P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

2.0 Material : Body – ABS

2.1 Color : Dark Blue


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com