

Technical Data Sheet

Doc No – Ab/TDS/Eqpt -129

Rev No : 00

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

SWIRLTOP- Ultra Flat Magnetic Stirrer

Valid for Catalog No. E11307

1. Description/Specification

1.1 Description: SWIRLTOP- Ultra Flat Magnetic Stirrer

1.1.1 This magnetic Stirrer is used in laboratory provides better mixing, Physical & Chemical analysis.

1.1.2 It comes with automatic reverse rotation for better mixing result.

1.1.3 Top plate is made from chemically resistant materials, slip proof & stable.

1.2 Design: Flat shaped Stirrer

1.3 Performance Specifications:

1.3.1 Speed Range : 15-1500 rpm.

1.3.2 Voltage : 220 ~ 240 V

1.3.3 Frequency : 50 Hz

1.3.4 Max. Stirrer Quantity : 0.8L

1.3.5 Max Load Capacity : 2 kg

1.3.6 Dimension of work place : 90x90 mm

1.3.7 Dimensions : 175x120x15 mm

1.3.8 Weight : 0.3 kg

1.3.9 Permissible ambient Temperature: 5-45°C

1.4.0 Permissible ambient : 80%

1.4 Material: Body – ABS Material

1.5 Color: Grey



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslalware.com
W : www.abdoslalware.com

Technical Data Sheet


1.6. Basic features:


1.6.1 **Material** : Made up of Highly pure Virgin ABS (Acrylonitrile Butadine Styrene)

1.6.2 **Dimension** : 175 X 120 X 15 mm

1.6.3 It contains special designed silicon feet base excellent for damping vibration.

1.6.4 Top plate is made from chemically resistant materials, slip proof & stable.


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttarakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata - 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com