

Technical Data Sheet

Doc No Ab/TDS/Eqpt-137

Rev No : 02

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

WAVEX – TUBE ROTATOR

Valid for Catlog No. E11270

1. Description/Specification

1.1 Description:

SWIRL Tube Rotator rotating Mixer gives gentle but effective mixing of biological samples to held in 1.5 ml to 50 ml tubes used in variety of application including immune precipitations, Prevention of blood coagulation and latex diagnostics.

1.2 Design :

Rotisserie Type

1.3 performance Specifications :

- 1.3.1 Voltage (AC) : 220~240 Volt
- 1.3.2 Frequency (Hz) : 50
- 1.3.3 Power (W) : 40
- 1.3.4 Angle Range : 0- 90
- 1.3.5 Motor type : DC Motor
- 1.3.6 Speed range (RPM) : 0 ~ 70
- 1.3.7 Speed display : LCD
- 1.3.8 Timer range (Min) : 1 – 1199
- 1.3.9 Timer display : LCD
- 1.4.0 Run type : Timer/Continuous operation
- 1.4.1 Overall dimension (mm) : 280 x 210 x 300 mm
- 1.4.2 Weight : 8.0 Kg
- 1.4.3 Permissible ambient Temp : 5- 45
- 1.4.4 Permissible Relative humidity : 80 %



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslaware.com
W : www.abdoslaware.com

Technical Data Sheet

1.4.5 Dual Dish accessories for large tubes

1.4.6 Protection class according to DIN EN60529 : IP21

1.4 Material : Body – ABS, Stainless Steel

1.5 Color : Dark Blue

2. Spare Accessories :


E11271 Disk accessory for 1.5 ml x 60 micro centrifuge tubes

E11272 Disk accessory for 50 ml x 8 Centrifuge tubes

E11273 Disk accessory for 15 ml x 16 tubes

E11274 Disk Support Rods used with disc accessory


Quality Assurance Incharge
Abdos Labtech Pvt. Ltd.



Shiv Ganga Industrial Estate,
Plot no. 1 Lakeshwari, Bhagwanpur,
Roorkee - 247667 Uttrakhand, India
An ISO 9001:2008 Company
A Division of Abdos

Regd. Office :
13A, Government Place (East)
Kolkata – 700069, India
CIN : U85100WB2008PTC123903

P : +91 (0) 1332 235277
E : info@abdoslabware.com
W : www.abdoslabware.com